

Canine/Feline Allergen-specific Antibody Rapid Test Kit (Colloidal gold)

Analysis of performance evaluation data

1. Determination method

The analytical performance indicators and test methods for the product are as follows:

1.1 Analytical performance

- 1.1.1 Physical inspection: The test cassette should be visually flat with firmly attached materials; components should be intact and in correct quantity; liquid components should show no obvious sedimentation or flocculation.
- 1.1.2 Blank limit: Repeat the detection of the buffer in the test kit, and all results should be negative.
- 1.1.3 Positive concordance rate: Each test item is evaluated using 1-3 positive quality control samples, and the results should match the characteristics of the positive quality control samples.
- 1.1.4 Negative concordance rate: A single negative quality control sample should yield negative results for all tests.
- 1.1.5 Precision (inter-laboratory variation): When testing the same batch of reagents with 1-3 positive quality control samples in three different laboratories, the results should match the characteristics of the positive quality control samples, and the color development should be uniform across repeated tests of the same quality control sample. If no positive quality control samples are available, a single negative quality control sample should be used for testing, and the results should match the characteristics of the negative quality control samples, with consistent color development across repeated tests of the same quality control sample.
- 1.1.6 Precision (inter-batch variation): Test different batches of reagents using 1-3 positive quality control samples. The results should match the characteristics of the positive quality control samples, and the color development should be uniform between repeated tests of the same quality control sample.
- 1.1.7 Validity period: After the expiration date of the kit or after storage at 37°C for 7 days, 14 days, and 20 days, the accuracy and precision should be tested using 1-3 positive quality control samples. The results should meet the requirements of the positive quality control samples, and the color development should be uniform between repeated tests of the same quality control sample.

2. Testing data

2.1 Blank line evaluation experiment

2.1.1 Experimental objective

This experiment was conducted to evaluate the blank limit of the company's produced Canine/Feline Allergen-specific Antibody Rapid Test Kit (Colloidal gold)

2.1.2 Experimental materials:

1. Feline inhalant allergen-specific antibody rapid test kit (Colloidal gold)
(Batch codes: KY25111105, KY25111315, KY25111515);
2. Canine inhalant allergen-specific antibody rapid test kit (Colloidal gold)
(Batch codes: KY25111106, KY25111316, KY25111516);
3. Feline food allergen-specific antibody rapid test kit (Colloidal gold)
(Batch codes: KY25111107, KY25111317, KY25111517);
4. Canine food allergen-specific antibody rapid test kit (Colloidal gold)
(Batch codes: KY25111108, KY25111318, KY25111518).

2.1.3 Experimental methods (n=12):

Perform the experiment according to the instructions. For the three batches of Feline/Canine allergen-specific antibody test kit (colloidal gold), use the included buffer in the kit for 3 days, repeating the test 4 times daily, resulting in a total of 12 data sets. Calculate the compliance rate of the blank limit.

2.1.4 Experimental Results

Kit name	Kit batch code	Blank line compliance rate
Feline inhalant allergen-specific antibody rapid test kit (Colloidal gold)	KY25111105	100%
	KY25111107	100%
	KY25111315	100%
Feline food allergen-specific antibody rapid test kit (Colloidal gold)	KY25111317	100%
	KY25111515	100%
	KY25111517	100%
Canine inhalant allergen-specific antibody rapid test kit (Colloidal gold)	KY25111106	100%
	KY25111108	100%
	KY25111316	100%
Canine food allergen-specific antibody rapid test kit (Colloidal gold)	KY25111318	100%
	KY25111516	100%
	KY25111518	100%

2.1.5 Experiment conclusion :

The above experimental results indicate that the blank limit compliance rate of the Canine/Feline allergen-specific antibody test kit (Colloidal gold) meets the requirements.

2.2 Positive concordance rate evaluation test

Based on the following precision data, the results meet the criteria of the positive quality control samples, with a positive agreement rate of 100%.

2.3 Negative concordance rate assessment test

See the precision data below. The results are consistent with the properties of the negative quality control sample, with a negative compliance rate of 100%.

2.4 Precision evaluation experiment

2.4.1 Experimental objective

This experiment was conducted to evaluate the precision of the colloidal gold method for detecting allergen-specific antibodies in feline and canine species.

2.4.2 Experimental materials

1. Feline inhalant allergen-specific antibody test kit (Batch No.: KY25111105, KY25111315, KY25111515);
2. Canine inhalant allergen-specific antibody test kit (Batch No.: KY25111106, KY25111316, KY25111516);
3. Canine inhalant allergen-specific antibody test kit (Batch No.: KY25111107, KY25111317, KY25111517);
4. Canine food allergen-specific antibody test kit (Batch No.: KY25111108, KY25111318, KY25111518).
5. Quality control materials

Quality control sample number	Quality control attributes	Properties of quality control materials
Cat 5	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +
Cat 7	Negative quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> - Soybean - Egg white - Beef - Shrimp - Milk -
Cat 8	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +
Cat 9	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +
Cat 10	Positive quality control	Mugwort - Candida albicans -

		Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +
Cat 13	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum + Cockroach - <i>Dermatophgoides farinae</i> +
Cat 24	Positive quality control	Mugwort - Candida albicans,- Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +

Quality control sample number	Quality control attributes	Properties of quality control materials
Dog 3	Positive quality control	Mugwort - Candida albicans + Penicillium chrysogenum + Cockroach - <i>Dermatophgoides farinae</i> +
Dog 60	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +
Dog 68	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum + Cockroach - <i>Dermatophgoides farinae</i> +
Dog 81	Negative quality control	Soybean - Egg white - Beef - Shrimp - Milk -
Dog 102	Negative quality control	Soybean - Egg white - Beef - Shrimp - Milk -
Dog 106	Positive quality control	Mugwort - Candida albicans + Penicillium chrysogenum +

		Cockroach - <i>Dermatophgoides farinae</i> +
Dog 113	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> +
Dog 147	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum +++ Cockroach - <i>Dermatophgoides farinae</i> +

2.4.3 Experimental methodology

Inter-laboratory variation assessment method: Use the Feline/Canine allergen-specific antibody test kit (Colloidal gold) to conduct 1-3 repeated tests on each positive sample (n=1-3). If no positive samples are obtained, one negative sample is used. This experiment utilized a batch of kits, with 5-day repeated testing across 3 laboratories, each test performed 5 times. A total of 75 replicate data points were obtained from each control sample.

Inter-batch variance assessment method: Repeat the detection of 1-3 positive specimens (n=1-3) using three different batch codes of feline/dog allergen-specific antibody test kits (colloidal gold method). If no positive samples are obtained, use one negative sample. Strictly follow the instructions to calculate inter-laboratory variance and inter-batch variance.

2.4.4 Experimental Results

2.4.4.1 Intertrial variation (n=3–9)

The specific experimental results are shown in the table below:

The precision (inter-laboratory variation) results of the **Feline inhalant** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111105 are as follows:

Test days	Positive quality control sample 1 (Cat 8)			Positive quality control sample 2 (Cat 10)			Positive quality control sample 3 (Cat 24)		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 2	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+

Note: 'T5+' indicates 'positive for *Dermatophgoides farinae* in the fifth test item of this panel,' while 'T5++' indicates 'positive for *Dermatophgoides farinae* in the fifth test item of this panel.'

The precision (inter-laboratory variation) results of the **Feline Food** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111107:

Test days	Negative quality control sample (Cat 7)		
	Lab 1	Lab 2	Lab 3
Day 1	-	-	-
Day 2	-	-	-
Day 3	-	-	-
Day 4	-	-	-
Day 5	-	-	-

Note: The symbol "-" indicates "the test result of this panel is negative", the same applies below.

The precision (inter-laboratory variation-Day 1) results of the **Canine inhalant** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111106 are as follows:

Allergen test items	Positive quality control sample 1 (Dog 106)			Positive quality control sample 2 (Dog 3)			Positive quality control sample 3 (Dog 68)		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	+	+	+	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	+	+	+
Cockroach	-	-	-	-	-	-	-	-	-
Dermatophgoides farinae	+	+	+	+	+	+	+	+	+

Note: The "+" indicates "positive for the allergen in this test panel," and the "-" indicates "negative" for the allergen in this test panel," with the same notation applied below.

The precision (inter-laboratory variation-Day 2) results of the **Canine inhalant** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111106 are as follows:

Allergen test items	Positive quality control sample 1 (Dog 106)			Positive quality control sample 2 (Dog 3)			Positive quality control sample 3 (Dog 68)		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	+	+	+	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	+	+	+
Cockroach	-	-	-	-	-	-	-	-	-
Dermatophgoides farinae	+	+	+	+	+	+	+	+	+

The precision (inter-laboratory variation-Day 3) results of the **Canine inhalant** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111106 are as follows:

Allergen test items	Positive quality control sample 1 (Dog 106)	Positive quality control sample 2 (Dog 3)	Positive quality control sample 3 (Dog 68)
---------------------	--	--	---

	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	+	+	+	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	+	+	+
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	+	+	+	+	+	+

The precision (inter-laboratory variation-**Day 4**) results of the **Canine inhalant** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111106 are as follows:

Allergen test items	Positive quality control sample 1 (Dog 106)			Positive quality control sample 2 (Dog 3)			Positive quality control sample 3 (Dog 68)		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	+	+	+	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	+	+	+
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	+	+	+	+	+	+

The precision (inter-laboratory variation-**Day 5**) results of the **Canine inhalant** allergen-specific antibody test kit (Colloidal gold) Batch No. KY25111106 are as follows:

Allergen test items	Positive quality control sample 1 (Dog 106)			Positive quality control sample 2 (Dog 3)			Positive quality control sample 3 (Dog 68)		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	+	+	+	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	+	+	+
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	+	+	+	+	+	+

The precision (inter-laboratory variation) results of the **Canine food** allergen-specific antibody test kit (Colloidal gold) Batch No. (KY25111108) are as follows:

Allergen test items	Negative quality control (Dog 81)														
	Day 1			Day 2			Day 3			Day 4			Day 5		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Soybean	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Egg white	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Beef	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Shrimp	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Milk	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

The experimental results demonstrate that Feline inhalant/food allergen-specific antibody test kit after excluding visual interpretation discrepancies and laboratory temperature effects, shows consistent results across different laboratories with matching quality control standards, and uniform color development in repeated tests using the same control sample. Similarly, the Canine inhalant/food allergen-specific antibody test kit (colloidal gold) also exhibits comparable performance, with consistent results across laboratories and uniform color development in repeated tests. These findings confirm that the inter-laboratory variability meets the required standards.

The experimental results show that both the Feline and Canine inhalant/food allergen-specific antibody test kits (colloidal gold) perform consistently across different laboratories. When visual interpretation discrepancies and laboratory temperature effects are excluded, both kits demonstrate matching quality control standards and uniform color development in repeated tests using the same control sample. These findings confirm that the inter-laboratory variability meets the required standards.

2.4.4.2 Inter-laboratory variation (n=30) ;

The specific experimental results are shown in the table below:

The precision (inter-batch difference) results of the **Feline inhalant** allergen-specific antibody test kit (Colloidal gold) are as follows:

Allergen test items	Positive quality control sample 1 (Cat 5)			PositiveQuality Control Sample 2 (Cat 9)			Positive quality control sample 3 (Cat 13)		
	KY25111108	KY25111315	KY25111515	KY25111108	KY25111315	KY25111515	KY25111108	KY25111315	KY25111515
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	-	-	-	-	-	-	-	-	-
Penicillium chrysogenum	-	-	-	-	-	-	+	+	+
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	+	+	+	+	+	+

The precision (inter-batch difference) results of the **Feline food** allergen-specific antibody test kit (Colloidal gold) are as follows:

Allergen test items	Negative quality control (Cat 7)		
	KY25111107	KY25111317	KY25111517
Soybean	-	-	-
Egg white	-	-	-
Beef	-	-	-
Shrimp	-	-	-
Milk	-	-	-

The precision (inter-batch difference) results of the **Canine inhalant** allergen-specific antibody test kit (Colloidal gold) are as follows:

Allergen test items	Positive quality control sample 1 (Dog 60)			Positive quality control sample 2 (Dog 113)			Positive quality control sample 3 (Dog 147)		
	KY25111106	KY25111316	KY25111516	KY25111106	KY25111316	KY25111516	KY25111106	KY25111316	KY25111516
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	-	-	-	+	+	+	-	-	-
Penicillium chrysogenum	-	-	-	+	+	+	+++	+++	+++
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	+	+	+	+	+	+

Note: The notation "+++" indicates "strongly positive for the allergens in this panel", the same applies hereafter.

The precision (inter-batch difference) results of the **Canine food** allergen-specific antibody test kit (Colloidal gold) are as follows:

Allergen test items	Negative quality control (dog 102)		
	KY25111108	KY25111318	KY25111518
Soybean	-	-	-
Egg white	-	-	-
Beef	-	-	-
Shrimp	-	-	-
Milk	-	-	-

The experimental results demonstrate that the colloidal gold-based test kit for feline inhalant/food allergen-specific antibody detection shows consistent results across batches, with uniform color development in repeated tests using the same control sample. Similarly, the canine version of this kit demonstrates batch-to-batch consistency and uniform color development in repeated tests with the same control sample. These findings confirm that the test kits meet the required standards without inter-batch variability.

2.5 Validity evaluation test

2.5.1 Experimental objective

This experiment was conducted to evaluate the shelf life of the company's Colloidal Gold Assay Kit for Specific Antibody Detection of Feline/Dog Allergens.

2.5.2 Experimental materials:

1. Feline inhalant allergen-specific antibody test kit (Colloidal gold) (Batch codes: KY25111105, KY25111315, KY25111515);
2. Feline food allergen-specific antibody test kit (Colloidal gold) (Batch codes: KY25111107, KY25111317, KY25111517);
3. Canine inhalant allergen-specific antibody test kit (Colloidal gold) (Batch codes.: KY25111106, KY25111316, KY25111516);
4. Canine food allergen-specific antibody test kit (Colloidal gold) (Batch codes: KY25111108, KY25111318, KY25111518);
5. Quality control materials

Quality control sample number	Quality control attributes	Properties of quality control materials
Cat 4	Negative quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> - Soybean - Egg white - Beef - Shrimp - Milk -
Cat 8	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> ++
Cat 10	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> ++
Cat 24	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> ++

Quality control sample number	Quality control attributes	Properties of quality control materials
Dog 14	Positive quality control	Mugwort - Candida albicans + Penicillium chrysogenum + Cockroach - <i>Dermatophgoides farinae</i> +++
Dog 27	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum + Cockroach - <i>Dermatophgoides farinae</i> +
Dog 99	Positive quality control	Mugwort - Candida albicans - Penicillium chrysogenum - Cockroach - <i>Dermatophgoides farinae</i> ++
Dog 166	Positive quality control	Soybean - Egg white Beef - Shrimp - Milk +

2.5.3 Experimental methods (n=12):

Perform experiments according to the instructions, test three batches of reagent kits, and store the kits at 37°C for 7 days, 14 days, and 20 days respectively. Conduct tests on the same day, with three replicates at each temperature. A total of three batches were tested, yielding nine data sets, to evaluate the shelf life of the test kits.

The specific experimental data results are as follows:

Feline inhalant allergen-specific antibody test kit (Colloidal gold); Batch No.:KY25111105; Test results as follows:

Allergen test items	Positive quality control sample 1 (Cat 8)			Positive quality control sample 2 (Cat 10)			Positive quality control sample 3 (Ca t24)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	-	-	-	-	-	-	-	-	-
Penicillium chrysogenum	-	-	-	-	-	-	-	-	-
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	++	++	++	+	+	+

Note: The notation "++" indicates "positive for the allergen in this panel", the same applies hereafter.

Feline inhalant allergen-specific antibody test kit (Colloidal gold); Batch No.: KY25111315; Test results as follows:

Allergen test items	Positive quality control sample 1 (Cat 8)			Positive quality control sample 2 (Cat 10)			Positive quality control sample 3 (Ca t24)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	-	-	-	-	-	-	-	-	-
Penicillium chrysogenum	-	-	-	-	-	-	-	-	-
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	++	++	++	+	+	+

Feline inhalant allergen-specific antibody test kit (Colloidal gold); Batch No.: KY25111515; Test results as follows:

Allergen test items	Positive quality control sample 1 (Cat 8)			Positive quality control sample 2 (Cat 10)			Positive quality control sample 3 (Ca t24)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	-	-	-	-	-	-	-	-	-
Penicillium chrysogenum	-	-	-	-	-	-	-	-	-
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+	+	+	++	++	++	+	+	+

Feline food allergen-specific antibody test kit (Colloidal gold); Batch No.:KY25111107 , KY25111317 , KY25111517 ; Test results as follows:

Allergen test items	KY25111107 Negative control (Cat 4)			KY25111317 Negative control (Cat 4)			KY25111517 Negative control (Cat 4)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Soybean	-	-	-	-	-	-	-	-	-
Egg white	-	-	-	-	-	-	-	-	-
Beef	-	-	-	-	-	-	-	-	-
Shrimp	-	-	-	-	-	-	-	-	-
Milk	-	-	-	-	-	-	-	-	-

Canine inhalant allergen-specific antibody test kit (Colloidal gold); Batch No.: KY25111106; Test results as follows:

Allergen test items	Positive quality control sample 1 (Dog 14)			Positive quality control sample 2 (Dog 27)			Positive quality control sample 3 (Dog 99)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	-	-	-	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	-	-	-
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+++	+++	+++	+	+	+	++	++	++

Canine inhalant allergen-specific antibody test kit (Colloidal gold); Batch No.: KY25111316; Test results as follows::

Allergen test items	Positive quality control sample 1 (Dog 14)			Positive quality control sample 2 (Dog 27)			Positive quality control sample 3 (Dog 99)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	-	-	-	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	-	-	-
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+++	+++	+++	+	+	+	++	++	++

Canine inhalant allergen-specific antibody test kit (Colloidal gold); Batch

No.:KY25111516; Test results as follows:

Allergen testing items	Positive quality control sample 1 (Dog 14)			Positive quality control sample 2 (Dog 27)			Positive quality control Sample 3 (Dog 99)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Mugwort	-	-	-	-	-	-	-	-	-
Candida albicans	+	+	+	-	-	-	-	-	-
Penicillium chrysogenum	+	+	+	+	+	+	-	-	-
Cockroach	-	-	-	-	-	-	-	-	-
<i>Dermatophgoides farinae</i>	+++	+++	+++	+	+	+	++	++	++

Canine food allergen-specific antibody test kit (Colloidal gold); Batch No.: KY25111108, KY25111318, KY25111518 ; Test results as follows:

Allergen test items	KY25111108 Positive control sample (Dog 166)			KY25111318 Positive control sample (Dog 166)			KY25111518 Positive control sample (Canine 166)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
Soybean	-	-	-	-	-	-	-	-	-
Egg white	-	-	-	-	-	-	-	-	-
Beef	-	-	-	-	-	-	-	-	-
Shrimp	-	-	-	-	-	-	-	-	-
Milk	+	+	+	+	+	+	+	+	+

The experimental results demonstrate that the test outcomes meet the requirements for both positive and negative quality control samples, with consistent color development across repeated tests of the same control sample. In accordance with the product technical specifications, the kit must be stored at 37°C for 20 days from the date of production (with a shelf life of 24 months). Consequently, The following kits have a shelf life of 24 months: Feline inhalant allergen-specific antibody test kit (Colloidal gold) , Feline food allergen-specific antibody test kit (Colloidal gold) , Canine inhalant allergen-specific antibody test kit (Colloidal gold), Canine food allergen-specific antibody test kit (Colloidal gold) .