

Analysis of Performance Evaluation Data for Canine/Feline Allergen-Specific Antibody Rapid Test Kit

1. Protocol Definition

The analytical performance indicators and test methods for the product are as follows:

1.1 Analytical performance Indicators

- 1.1.1 **Physical inspection:** The test cassette should be visually flat with firmly attached materials; components should be intact and in correct quantity; liquid components should show no obvious sedimentation or flocculation.
- 1.1.2 **Limit of Blank (LoB):** Repeat the buffer in the test kit, and all results should be negative.
- 1.1.3 **Positive Percent Agreement (PPA):** Each test item is evaluated using 1-3 positive controls, and the results should match the property of the positive controls.
- 1.1.4 **Negative Percent Agreement (NPA):** A single negative control should yield negative results for all tests.
- 1.1.5 **Precision (inter-laboratory variation):** When testing the same batch of kits with 1-3 positive controls in three different laboratories, the results should match the properties of the positive controls, and the color development should be uniform across repeated tests of the same quality control sample. If no positive controls are available, a single negative control should be used for testing, and the results should match the properties of the negative controls, with consistent color development across repeated tests of the same quality control sample.
- 1.1.6 **Precision (inter-batch variation):** Test different batches of kit using 1-3 positive controls. The results should match the characteristics of the positive controls, and the color development should be uniform between repeated tests of the same quality control sample.
- 1.1.7 **Shelf-Life (Accelerated Stability):** After the expiration date of the kit or after storage at 37°C for 7 days, 14 days, and 20 days, the accuracy and precision should be tested using 1-3 positive controls. The results should meet the requirements of the positive controls, and the color development should be uniform between repeated tests of the same quality control sample.

2. Testing data

2.1 Limit of blank (LoB) evaluation experiment

2.1.1 Experimental objective

This experiment was conducted to evaluate the limit of blank (LoB) of the company's canine/feline allergen-specific antibody rapid test kits (multiple panels and multiple lots).

2.1.2 Experimental materials:

Product name	Catalogue No.	Kit batch codes
Canine allergen-specific	G-DOG-ISC01	KY25111106, KY25111316, KY25111516
	G-DOG-ISC02	KY25110602, KY25111202, KY25111402

antibody rapid test kit	G-DOG-ISC03	KY25110606, KY25111206, KY25111406
	G-DOG-ISC04	KY25110608, KY25111208, KY25111408
	G-DOG-ISC05	KY25110502, KY25111210, KY25111410
	G-DOG-ISC06	KY25110610, KY25111212, KY25111412
	G-DOG-FSC01	KY25111108, KY25111318, KY25111518
	G-DOG-FSC02	KY25110604, KY25111204, KY25111404
	G-DOG-FSC03	KY25110606, KY25111206, KY25111406
	G-DOG-FSC04	KY25111002, KY25111216, KY25111416
	G-DOG-FSC05	KY25110702, KY25111218, KY25111418
	G-DOG-FSC06	KY25110704, KY25111302, KY25111502
	G-DOG-FSC07	KY25110706, KY25111304, KY25111504
	G-DOG-FSC08	KY25111004, KY25111306, KY25111506
	G-DOG-MIX01	KY25111006, KY25111308, KY25111508
	G-DOG-MIX02	KY25111008, KY25111310, KY25111510
	G-DOG-MIX03	KY25111102, KY25111312, KY25111512
	G-DOG-MIX04	KY25111104, KY25111314, KY25111514

Product name	Catalogue No.	Kit batch codes
Feline allergen-specific antibody rapid test kit	G-CAT-ISC01	KY25111105, KY25111315, KY25111515
	G-CAT-ISC02	KY25110601, KY25111201, KY25111401
	G-CAT-ISC02	KY25110601, KY25111201, KY25111401
	G-CAT-ISC03	KY25110605, KY25111205, KY25111405
	G-CAT-ISC04	KY25110607, KY25111207, KY25111407
	G-CAT-ISC05	KY25110501, KY25111209, KY25111409
	G-CAT-ISC06	KY25110609, KY25111211, KY25111411
	G-CAT-FSC01	KY25111107, KY25111317, KY25111517
	G-CAT-FSC02	KY25110603, KY25111203, KY25111403
	G-CAT-FSC03	KY25110707, KY25111213, KY25111413
	G-CAT-FSC04	KY25111001, KY25111215, KY25111415
	G-CAT-FSC05	KY25110701, KY25111217, KY25111417
	G-CAT-FSC06	KY25110703, KY25111301, KY25111501
	G-CAT-FSC07	KY25110705, KY25111303, KY25111503
	G-CAT-FSC08	KY25111003, KY25111305, KY25111505
	G-CAT-MIX01	KY25111005, KY25111307, KY25111507
	G-CAT-MIX02	KY25111007, KY25111309, KY25111509
	G-CAT-MIX03	KY25111101, KY25111311, KY25111511
	G-CAT-MIX04	KY25111103, KY25111313, KY25111513

2.1.3 Experimental methods (n=12):

Perform the experiment according to the instructions. For the three batches of Feline/Canine allergen-specific antibody rapid test kit, use the included buffer in the kit for 3 days, repeating the test 4 times daily, resulting in a total of 12 data sets. Calculate the compliance rate of the Limit of Blank (LoB).

2.1.4 Experimental Results

Product name	Catalogue No.	Kit batch codes	LoB compliance rate
Canine allergen-specific antibody rapid test kit	G-DOG-ISC01	KY25111106	100%
		KY25111516	100%
		KY25111316	100%
	G-DOG-ISC02	KY25110602	100%
		KY25111202	100%
		KY25111402	100%
	G-DOG-ISC03	KY25110606	100%
		KY25111206	100%
		KY25111406	100%
	G-DOG-ISC04	KY25110608	100%
		KY25111208	100%
		KY25111408	100%
	G-DOG-ISC05	KY25110502	100%
		KY25111210	100%
		KY25111410	100%
	G-DOG-ISC06	KY25110610	100%
		KY25111212	100%
		KY25111412	100%
	G-DOG-FSC01	KY25111318	100%
		KY25111108	100%
		KY25111518	100%
	G-DOG-FSC02	KY25110604	100%
		KY25111204	100%
		KY25111404	100%
	G-DOG-FSC03	KY25110606	100%
		KY25111206	100%
		KY25111406	100%
	G-DOG-FSC04	KY25111002	100%
		KY25111216	100%
		KY25111416	100%
	G-DOG-FSC05	KY25110702	100%
		KY25111218	100%
		KY25111418	100%
	G-DOG-FSC06	KY25110704	100%
		KY25111302	100%
		KY25111502	100%
	G-DOG-FSC07	KY25110706	100%
		KY25111304	100%
		KY25111504	100%
	G-DOG-FSC08	KY25111004	100%

		KY25111306	100%
		KY25111506	100%
	G-DOG-MIX01	KY25111006	100%
		KY25111308	100%
		KY25111508	100%
	G-DOG-MIX02	KY25111008	100%
		KY25111310	100%
		KY25111510	100%
	G-DOG-MIX03	KY25111102	100%
		KY25111312	100%
		KY25111512	100%
	G-DOG-MIX04	KY25111104	100%
		KY25111314	100%
		KY25111514	100%

Product name	Catalogue No.	Kit batch codes	LoB compliance rate
Feline allergen-specific antibody rapid test kit	G-CAT-ISC01	KY25111105	100%
		KY25111515	100%
		KY25111315	100%
	G-CAT-ISC02	KY25110601	100%
		KY25111201	100%
		KY25111401	100%
	G-CAT-ISC03	KY25110605	100%
		KY25111205	100%
		KY25111405	100%
	G-CAT-ISC04	KY25110607	100%
		KY25111207	100%
		KY25111407	100%
	G-CAT-ISC05	KY25110501	100%
		KY25111209	100%
		KY25111409	100%
	G-CAT-ISC06	KY25110609	100%
		KY25111211	100%
		KY25111411	100%
	G-CAT-FSC01	KY25111317	100%
		KY25111107	100%
		KY25111517	100%
	G-CAT-FSC02	KY25110603	100%
		KY25111203	100%
		KY25111403	100%
	G-CAT-FSC03	KY25110707	100%
		KY25111213	100%
		KY25111413	100%

	G-CAT-FSC04	KY25111001	100%
		KY25111215	100%
		KY25111415	100%
	G-CAT-FSC05	KY25110701	100%
		KY25111217	100%
		KY25111417	100%
	G-CAT-FSC06	KY25110703	100%
		KY25111301	100%
		KY25111501	100%
	G-CAT-FSC07	KY25110705	100%
		KY25111303	100%
		KY25111503	100%
	G-CAT-FSC08	KY25111003	100%
		KY25111305	100%
		KY25111505	100%
	G-CAT-MIX01	KY25111005	100%
		KY25111307	100%
		KY25111507	100%
	G-CAT-MIX02	KY25111007	100%
		KY25111309	100%
		KY25111509	100%
	G-CAT-MIX03	KY25111101	100%
		KY25111311	100%
		KY25111511	100%
	G-CAT-MIX04	KY25111103	100%
		KY25111313	100%
		KY25111513	100%

2.1.5 Experiment conclusion :

The above experimental results indicate that the blank limit compliance rate of the Canine/Feline allergen-specific antibody rapid test kit meets the requirements.

2.2 Positive percent agreement evaluation test

Based on the following precision data, the results meet the criteria of the positive controls, with a positive agreement rate of 100%.

2.3 Negative percent agreement assessment test

See the precision data below. The results are consistent with the properties of the negative control, with a negative compliance rate of 100%.

2.4 Precision evaluation experiment

2.4.1 Experimental objective

This experiment was conducted to evaluate the precision of the canine/Feline allergen-specific antibody rapid test kit

2.4.2 Experimental materials

Product name	Catalogue No.	Kit batch codes
Canine allergen-specific antibody rapid test kit	G-DOG-ISC01	KY25111106, KY25111316, KY25111516
	G-DOG-ISC02	KY25110602, KY25111202, KY25111402
	G-DOG-ISC03	KY25110606, KY25111206, KY25111406
	G-DOG-ISC04	KY25110608, KY25111208, KY25111408
	G-DOG-ISC05	KY25110502, KY25111210, KY25111410
	G-DOG-ISC06	KY25110610, KY25111212, KY25111412
	G-DOG-FSC01	KY25111108, KY25111318, KY25111518
	G-DOG-FSC02	KY25110604, KY25111204, KY25111404
	G-DOG-FSC03	KY25110606, KY25111206, KY25111406
	G-DOG-FSC04	KY25111002, KY25111216, KY25111416
	G-DOG-FSC05	KY25110702, KY25111218, KY25111418
	G-DOG-FSC06	KY25110704, KY25111302, KY25111502
	G-DOG-FSC07	KY25110706, KY25111304, KY25111504
	G-DOG-FSC08	KY25111004, KY25111306, KY25111506
	G-DOG-MIX01	KY25111006, KY25111308, KY25111508
	G-DOG-MIX02	KY25111008, KY25111310, KY25111510
	G-DOG-MIX03	KY25111102, KY25111312, KY25111512
	G-DOG-MIX04	KY25111104, KY25111314, KY25111514

Product name	Catalogue No.	Kit batch codes
Feline allergen-specific antibody rapid test kit	G-CAT-ISC02	KY25110601, KY25111201, KY25111401
	G-CAT-ISC03	KY25110605, KY25111205, KY25111405
	G-CAT-ISC04	KY25110607, KY25111207, KY25111407
	G-CAT-ISC05	KY25110501, KY25111209, KY25111409

	G-CAT-ISC06	KY25110609, KY25111211, KY25111411
	G-CAT-FSC02	KY25110603, KY25111203, KY25111403
	G-CAT-FSC03	KY25110707, KY25111213, KY25111413
	G-CAT-FSC04	KY25111001, KY25111215, KY25111415
	G-CAT-FSC05	KY25110701, KY25111217, KY25111417
	G-CAT-FSC06	KY25110703, KY25111301, KY25111501
	G-CAT-FSC07	KY25110705, KY25111303, KY25111503
	G-CAT-FSC08	KY25111003, KY25111305, KY25111505
	G-CAT-MIX01	KY25111005, KY25111307, KY25111507
	G-CAT-MIX02	KY25111007, KY25111309, KY25111509
	G-CAT-MIX03	KY25111101, KY25111311, KY25111511
	G-CAT-MIX04	KY25111103, KY25111313, KY25111513

1. Quality control materials

Control ID	Quality control attributes	Properties of quality control materials
Cat 1	Negative Control	White ash(-), Mulberry(-), Cottonwood(-), Willow(-), Sycamore(-), Cocklebur(-), Nettle(-), Japanese hop(-), Common ragweed(-), Mugwort(-), Wheat(-), Rice(-), Corn(-), Soybean(-), Sesame(-), Peanut(-), Cashew(-), Almond(-), Pistachio(-), Walnut(-), Apple(-), Pear(-), Pineapple(-), Watermelon(-), Kiwi(-)
Cat 5	Positive Control	<i>Candida albicans</i> (-), <i>Penicillium chrysogenum</i> (-), <i>Cockroach</i> (-), <i>Alternaria alternata</i> (-), <i>Malassezia</i> (-), <i>Tyrophagus putrescentiae</i> (-), <i>Blomia tropicalis</i> (-), <i>D. farinae</i> (+), Flea(-), Mosquito(-), Chicken feather(-), Duck feather(-), Goose feather(-), White ash(-), Mulberry(-), Cottonwood(-), Willow(-), Sycamore(-), Cocklebur(-), Nettle(-), Japanese hop(-), Common ragweed(-), <i>Mugwort</i> (-), Pumpkin(-), Potato(-), Carrot(-), Cauliflower(-), Tomato(-), Apple(-), Pear(-), Pineapple(-), Watermelon(-), Kiwi(-), Wheat(-), Salmon(-), Duck(-), Shrimp(-), Turkey(-), <i>Malassezia</i> (-)Willow(-)
Cat 7	Negative Control	<i>Soybean</i> (-), <i>Egg white</i> (-), <i>Beef -Shrimp</i> (-), <i>Milk</i> (-), <i>Candida albicans</i> (-), <i>Penicillium chrysogenum</i> (-), <i>Cockroach</i> (-), <i>D. farinae</i> (-), Flea(-), Mosquito(-), Chicken feather(-), Duck feather(-), Goose feather(-), White ash(-), Mulberry(-),

		Cottonwood(-), Willow(-), Sycamore(-), Cocklebur(-), Nettle(-), Japanese hop(-), Common ragweed(-), Mugwort(-) , Soybean(-), Chicken(-), Duck(-), Sardine(-), Salmon(-), Lamb(-), Pork(-), Turkey(-), Shrimp(-), Scallop(-), Squid(-), Sardine(-), Wheat(-), Rice(-), Corn(-), Sesame(-), Peanut(-), Cashew(-), Almond(-), Pistachio(-), Walnut(-), Pumpkin(-), Potato(-), Carrot(-), Cauliflower(-), Tomato(-), Apple(-), Pear(-), Pineapple(-), Watermelon(-), Kiwi(-), Blomia tropicalis(-), Malassezia(-), Alternaria alternata(-)
Cat 8	Positive Control	Mugwort(-), Candida albicans(-), Penicillium chrysogenum(-), Cockroach(-), Tree grass pollen(-), Mold(-), Flea(-), Animal dander(-), D. farinae.(+) , Soybean(-), Sardine(-), Chicken(-), Flea(-)
Cat 9	Positive Control	Mugwort(-), Candida albicans(-), Penicillium chrysogenum(-), Cockroach(-), Alternaria alternata(-), Malassezia(-), Tyrophagus putrescentiae(-), Blomia tropicalis(-), D. farinae.(+) , Flea(-), Mosquito(-), Chicken feather(-), Duck feather(-), Goose feather(-), Pumpkin(-), Potato(-), Carrot(-), Cauliflower(-), Tomato(-), Wheat(-), Salmon(-), Duck(-), Shrimp(-), Turkey(-), Common ragweed(-), Willow(-)
Cat 10	Positive Control	Mugwort(-), Candida albicans(-), Penicillium chrysogenum(-), Cockroach(-), Tree grass pollen(-), Mold(-), Flea(-), Animal dander(-), D. farinae.(+) , Soybean(-), Sardine(-), Chicken(-)
Cat 13	Positive Control	Mugwort(-), Candida albicans(-), Penicillium chrysogenum(-), Cockroach(-), Alternaria alternata(-), Malassezia(-), Tyrophagus putrescentiae(-), Blomia tropicalis(-), D. farinae.(+) , Flea(-), Mosquito(-), Chicken feather(-), Duck feather(-), Goose feather(-), Pumpkin(-), Potato(-), Carrot(-), Cauliflower(-), Tomato(-), Wheat(-), Salmon(-), Duck(-), Shrimp(-), Turkey(-), Common ragweed(-), Willow(-)
Cat 14	Negative Control	White ash(-), Mulberry(-), Cottonwood(-), Willow(-), Sycamore(-), Cocklebur(-), Nettle(-), Japanese hop(-), Common ragweed(-), Mugwort(-), Lamb(-), Pork(-), Turkey(-), Duck(-), Chicken(-), Shrimp(-), Scallop(-), Squid(-), Sardine(-), Salmon(-), Wheat(-), Rice(-), Corn(-), Soybean(-), Sesame(-), Peanut(-), Cashew(-), Almond(-), Pistachio(-), Walnut(-), Apple(-), Pear(-), Pineapple(-), Watermelon(-), Kiwi(-)
Cat 17	Negative	Lamb(-), Pork(-), Turkey(-), Duck(-), Chicken(-), Shrimp(-), Scallop(-), Squid(-), Sardine(-), Salmon(-)

	Control	
Cat 20	Negative Control	Soybean(-), Chicken(-), Duck(-), Sardine(-), Salmon(-), Squid(-), Lamb(-), Alternaria alternata(-), Mugwort(-), Sycamore(-)
Cat 21	Negative Control	Soybean(-), Chicken(-), Duck(-), Sardine(-), Salmon(-), Squid(-), Lamb(-), Alternaria alternata(-), Mugwort(-), Sycamore(-)
Cat 24	Positive Control	Mugwort(-), Candida albicans(-), Penicillium chrysogenum(-), Cockroach(-), Tree grass pollen(-), Mold(-), Flea(-), Animal dander(-), D. farinae(+), Soybean(-), Sardine(-), Chicken(-)

Quality control sample #	Quality control attributes	Properties of quality control materials
Dog 2	Positive Control	White ash(-), Mulberry(-), Cottonwood(-), Willow(-), Sycamore(-), Cocklebur(-), Nettle(-), Japanese hop(+), Common ragweed(+), Mugwort(-), Shrimp(-), Scallop(-), Squid(-), Sardine(-), Salmon(-), Turkey(-), Malassezia(-)
Dog 3	Positive quality control	Mugwort (-), Candida albicans (+), Penicillium chrysogenum (+), Cockroach (-), D. farinae (+)
Dog 5	Negative Control	Wheat(-), Rice(-), Corn(-), Soybean(-), Sesame(-), Shrimp(-), Turkey(-), Malassezia(-), Common ragweed(-), Willow(-), Squid(-), Lamb(-), Alternaria alternata(-), Mugwort(-), Sycamore(-)
Dog 16	Negative Control	White ash(-), Mulberry(-), Cottonwood(-), Willow(-), Sycamore(-), Cocklebur(-), Nettle(-), Japanese hop(-), Common ragweed(-), Mugwort(-)
Dog 60	Positive Control	Mugwort (-), Candida albicans (+), Penicillium chrysogenum (+), Cockroach (-), Tree grass pollen(-), Mold(-), Flea(-), Animal dander(-), D. farinae(+), Alternaria alternata(-), Malassezia(-), Tyrophagus putrescentiae(-), Blomia tropicalis(-), Soybean(-), Sardine(+), Chicken(-)
Dog 61	Positive Control	Flea(-), Mosquito(-), Chicken feather(-), Duck feather(-), Goose feather(-), Pumpkin(-), Potato(-), Carrot(-), Cauliflower(-), Tomato(-)
Dog 68	Positive quality control	Mugwort (-), Candida albicans (-), Penicillium chrysogenum (+), Cockroach (-), D. farinae (+)
Dog 81	Negative Control	Egg white(-), Beef (-), Shrimp(-), Milk(-), Soybean(-), Chicken(-), Duck(-), Sardine(-), Salmon(-)
Dog 89	Negative Control	Lamb(-), Pork(-), Turkey(-), Duck(-), Chicken(-), Shrimp(-), Scallop(-), Squid(-), Sardine(-), Salmon(-)

Dog 102	Negative Control	Soybean(-), Egg white(-), Beef (-), Shrimp(-), Milk(-), Peanut(-), Cashew(-), Almond(-), Pistachio(-), Walnut(-), Wheat(-), Salmon(-), Duck(-), Mosquito(-), Blomia tropicalis(+)
Dog 105	Negative Control	Wheat(-), rice(-), Corn(-), Soybean(-), Sesame(-), Peanut(-), Cashew(-), Almond(-), Pistachio(-), Walnut(-)
Dog 106	Positive Control	Mugwort (-), Candida albicans (+), Penicillium chrysogenum (+), Cockroach (-), D. farinae (+)
Dog 113	Positive Control	Mugwort (-), Candida albicans (-), Penicillium chrysogenum (-), Cockroach (-), D. farinae (+)
Dog 130	Positive Control	Soybean(+), Sardine(+), Chicken(-), Flea(-), D. farinae(+)
Dog 139	Positive Control	Tree grass pollen(-), Mold(+), Flea(-), Animal dander(-), D. farinae(+), Alternaria alternata(-), Malassezia(-), Tyrophagus putrescentiae(-), Blomia tropicalis(-)
Dog 147	Positive quality control	Mugwort (-), Candida albicans (-), Penicillium chrysogenum (+++), Cockroach (-), D. farinae (+)
Dog 155	Negative Control	Soybean(-), Chicken(-), Duck(-), Sardine(-), Salmon(-), Lamb(-), Pork(-), Turkey(-)
Dog 158	Negative Control	Flea(-), Mosquito(-), Chicken feather(-), Duck feather(-), Goose feather(-), Pumpkin(-), Potato(-), Carrot(-), Cauliflower(-), Tomato(-), Apple(-), Pear(-), Pineapple(-), Watermelon(-), Kiwi(+), Squid(-), Lamb(-), Alternaria alternata(-), Mugwort(-), Sycamore(-)

2.4.3 Experimental methodology

Inter-laboratory variation assessment method: Use the Feline/Canine allergen-specific antibody test kit to conduct 1-3 repeated tests on each positive sample (n=1-3). If no positive samples are obtained, one negative sample is used. This experiment utilized a batch of kits, with 5-day repeated testing across 3 laboratories, each test performed 5 times. A total of 75 replicate data points were obtained from each control sample.

Inter-batch variance assessment method: Repeat the detection of 1-3 positive specimens (n=1-3) using three different batch codes of feline/dog allergen-specific antibody test kits (colloidal gold method). If no positive samples are obtained, use one negative sample. Strictly follow the instructions to calculate inter-laboratory variance and inter-batch variance.

2.4.4 Experimental Results

2.4.4.1 Intertrial variation (n=3–9)

The specific experimental results are shown in the table below:

1. The precision (inter-laboratory variation) results of the **Feline** allergen-specific antibody rapid test kit G-CAT-IS01 are as follows:

Test days	Positive Control (Cat 8)			Positive Control (Cat 10)			Positive Control (Cat 24)		
	KY25111105			KY25111105			KY25111105		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 2	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+

Note: 'T5+' indicates 'positive for *D. farinae* in the fifth test item of this panel,' while 'T5++' indicates 'positive for *D. farinae* in the fifth test item of this panel.'

2. The precision (inter-laboratory variation) results of the **Feline** allergen-specific antibody rapid test kit G-CAT-ISC02 are as follows:

Test days	Positive Control (Cat 8)			Positive Control (Cat 8)			Positive Control (Cat 8)		
	KY25110601			KY25111201			KY25111401		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 2	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+

Note: 'T5+' indicates 'positive for *D. farinae* in the fifth test item of this panel,' while all other test items are negative. The same applies below.

3. The precision (inter-laboratory variation) results of the **Feline** allergen-specific antibody rapid test kit G-CAT-ISC03 are as follows:

Test days	Positive Control (Cat 9)			Positive Control (Cat 9)			Positive Control (Cat 9)		
	KY25110605			KY25111205			KY25111405		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+

Day 2	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+

4. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-ISC04** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25110607			KY25111207			KY25111407		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

Note: The symbol "-" indicates "the test result of this panel is negative", the same applies below.

5. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-ISC05** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25110501			KY25111209			KY25111409		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

6. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-ISC06** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
-----------	--------------------------	--	--	--------------------------	--	--	--------------------------	--	--

	KY25110609			KY25111211			KY25111411		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

7. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-FSC01** are as follows:

Test days	Negative Control (Cat 7)														
	KY25111107														
	Lab 1					Lab 2					Lab 3				
	Soybean	Egg white	Beef	Shrimp	Milk	Soybean	Egg white	Beef	Shrimp	Milk	Soybean	Egg white	Beef	Shrimp	Milk
Day 1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

8. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-FSC02** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25110603			KY25111203			KY25111403		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-

Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

9. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-FSC03** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25110707			KY25111213			KY25111413		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

10. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-FSC04** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25111001			KY25111215			KY25111415		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

11. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit **G-CAT-FSC05** are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25110701			KY25111217			KY25111417		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3

Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

12. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC06 are as follows:

Test days	Positive Control sample (Cat 7)			Positive Control sample (Cat 7)			Positive Control sample (Cat 7)		
	KY25110703			KY25111301			KY25111501		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

13. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC07 are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25110705			KY25111303			KY25111503		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

14. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC08 are as follows:

Test days	Negative Control (Cat 7)		Negative Control (Cat 7)		Negative Control (Cat 7)	
-----------	--------------------------	--	--------------------------	--	--------------------------	--

	KY25111003			KY25111305			KY25111505		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

15. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX01 are as follows:

Test days	Positive Control (Cat 24)			Positive Control (Cat 24)			Positive Control (Cat 24)		
	KY25111005			KY25111307			KY25111507		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 2	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+

16. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX02 are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25111007			KY25111309			KY25111509		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

17. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX03 are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25111101			KY251111311			KY251111511		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

18. The precision (inter-laboratory variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX04 are as follows:

Test days	Negative Control (Cat 7)			Negative Control (Cat 7)			Negative Control (Cat 7)		
	KY25111103			KY251111313			KY251111513		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

19. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC01 are as follows:

Test days	Positive Control (Dog 106)			Positive Control (Dog 3)			Positive Control (Dog 68)		
	KY25111106			KY25111106			KY25111106		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T3+ T5+	T3+ T5+	T3+ T5+
Day 2	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T3+ T5+	T3+ T5+	T3+ T5+
Day 3	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T3+ T5+	T3+ T5+	T3+ T5+
Day 4	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T3+ T5+	T3+ T5+	T3+ T5+

Day 5	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T2+ T3+ T5+	T3+ T5+	T3+ T5+	T3+ T5+
-------	-------------	-------------	-------------	-------------	-------------	-------------	---------	---------	---------

Note: Catalogue No.: G-DOG-ISC01 's detectable allergens:

T1 Mugwort, T2 Candida albicans, T3 Penicillium chrysogenum, T4 Cockroach, T5 D. farinae.

20. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC02 are as follows:

Test days	Positive Control (Dog 139)			Positive Control (Dog 139)			Positive Control (Dog 139)		
	KY25110602			KY25111202			KY25111402		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+
Day 2	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+
Day 3	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+
Day 4	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+
Day 5	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+	T2+T5+

21. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC03 are as follows:

Test days	Positive Control (Dog 139)			Positive Control (Dog 139)			Positive Control (Dog 139)		
	KY25110606			KY25111206			KY25111406		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 2	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+

22. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC04 are as follows:

Test days	Negative Control (Dog 61)			Negative Control (Dog 61)			Negative Control (Dog 61)		
	KY25110608			KY25111208			KY25111408		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3

Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

23. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC05 are as follows:

Test days	Negative Control (Dog 16)			Negative Control (Dog 16)			Negative Control (Dog 16)		
	KY25110502			KY25111210			KY25111410		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

24. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC06 are as follows:

Test days	Negative Control (Dog 16)			Negative Control (Dog 16)			Negative Control (Dog 16)		
	KY25110610			KY25111212			KY25111412		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

25. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC01 are as follows:

Test days	Negative Control (Dog 81)														
	KY25111108														
	Lab 1					Lab 2					Lab 3				
	Soybean	Egg white	Beef	Shrimp	Milk	Soybean	Egg white	Beef	Shrimp	Milk	Soybean	Egg white	Beef	Shrimp	Milk
Day 1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

26. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC02 are as follows:

Test days	Negative Control (Dog 81)			Negative Control (Dog 81)			Negative Control (Dog 81)		
	KY25110604			KY25111204			KY25111404		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

27. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC03 are as follows:

Test days	Negative Control (Dog 89)			Negative Control (Dog 89)			Negative Control (Dog 89)		
	KY25110606			KY25111206			KY25111406		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-

Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

28. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC04 are as follows:

Test days	Negative Control (Dog 89)			Negative Control (Dog 89)			Negative Control (Dog 89)		
	KY25111002			KY25111216			KY25111416		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

29. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC05 are as follows:

Test days	Negative Control (Dog 105)			Negative Control (Dog 105)			Negative Control (Dog 105)		
	KY25110702			KY25111218			KY25111418		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

30. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC06 are as follows:

Test days	Negative Control (Dog 105)			Negative Control (Dog 105)			Negative Control (Dog 105)		
	KY25110704			KY25111302			KY25111502		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3

Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

31. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC07 are as follows:

Test days	Negative Control (Dog 61)			Negative Control (Dog 61)			Negative Control (Dog 61)		
	KY25110706			KY25111304			KY25111504		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+
Day 2	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+
Day 3	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+
Day 4	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+
Day 5	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+	T4+

32. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC08 are as follows:

Test days	Negative Control (Dog 158)			Negative Control (Dog 158)			Negative Control (Dog 158)		
	KY25111004			KY25111306			KY25111506		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

33. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX01 are as follows:

Test days	Positive Control (Dog 130)			Positive Control (Dog 130)			Positive Control (Dog 130)		
	KY25111006			KY25111308			KY25111508		

	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++
Day 2	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++
Day 3	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++
Day 4	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++
Day 5	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++	T1+T2+T5++

Note: "T2+" means "the second test item of the panel - sardine allergen is positive", and "T5++" means "the fifth test item of the panel - *D. farinae* allergen is positive".

34. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX02 are as follows:

Test days	Positive Control (Dog 102)			Positive Control (Dog 102)			Positive Control (Dog 102)		
	KY25111008			KY25111310			KY25111510		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 2	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 3	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 4	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+
Day 5	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+	T5+

35. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX03 are as follows:

Test days	Negative Control (Dog 5)			Negative Control (Dog 5)			Negative Control (Dog 5)		
	KY25111102			KY25111312			KY25111512		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

36. The precision (inter-laboratory variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX04 are as follows:

Test days	Negative Control (Dog 5)			Negative Control (Dog 5)			Negative Control (Dog 5)		
	KY25111104			KY25111314			KY25111514		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
Day 1	-	-	-	-	-	-	-	-	-
Day 2	-	-	-	-	-	-	-	-	-
Day 3	-	-	-	-	-	-	-	-	-
Day 4	-	-	-	-	-	-	-	-	-
Day 5	-	-	-	-	-	-	-	-	-

The experimental results demonstrate that feline allergen-specific antibody rapid test kit (after excluding visual interpretation discrepancies and laboratory temperature effects, shows consistent results across different laboratories with matching quality control standards, and uniform color development in repeated tests using the same control sample. Similarly, the canine allergen-specific antibody rapid test kit also exhibits comparable performance, with consistent results across laboratories and uniform color development in repeated tests. These findings confirm that the inter-laboratory variability meets the required standards.

2.4.4.2 Inter-batch variation (n=30) ;

The specific experimental results are shown in the table below:

1. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-ISC01 are as follows:

Kit batch code	Positive Control (Cat 5)	Positive Control (Cat 9)	Positive Control (Cat 13)
KY25111108	T5+	T5+	T3+ T5+
KY25111315	T5+	T5+	T3+ T5+
KY25111515	T5+	T5+	T3+ T5+

Note: T3: *Penicillium chrysogenum* T5: *D. farinae*

2. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-ISC02 are as follows:

Kit batch code	Positive Control (Cat24)	Positive Control (Cat8)	Positive Control (Cat10)
KY25110601	T5+	T5+	T5+
KY25111201	T5+	T5+	T5+
KY25111401	T5+	T5+	T5+

3. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-ISC03 are as follows:

Kit batch code	Positive Control (Cat5)	Positive Control (Cat9)	Positive Control (Cat13)
KY25110605	T5+	T5+	T5+
KY25111205	T5+	T5+	T5+
KY25111405	T5+	T5+	T5+

4. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-ISC04 are as follows:

Kit batch code	Negative Control (Cat 5)	Negative Control (Cat 9)	Negative Control (Cat 13)
KY25110607	-	-	-
KY25111207	-	-	-
KY25111407	-	-	-

5. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-ISC05 are as follows:

Kit batch code	Negative Control (Cat5)	Negative Control (Cat 1)	Negative Control (Cat 14)
KY25110501	-	-	-
KY25111209	-	-	-
KY25111409	-	-	-

6. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-ISC06 are as follows:

Kit batch code	Negative Control (Cat 5)	Negative Control (Cat 1)	Negative Control (Cat1 4)
KY25110609	-	-	-
KY25111211	-	-	-
KY25111411	-	-	-

7. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC01 are as follows:

3. Kit batch code	4. Negative Control (Cat 7)	5. Negative Control (Cat 7)	6. Negative Control (Cat 7)
7. KY25111107	8. -	9. -	10. -
11. KY25111317	12. -	13. -	14. -
15. KY25111517	16. -	17. -	18. -

Note : G-CAT-FSC01 T1 Soybean T2 Egg white T3 Beef, T4 Shrimp T5 Milk

8. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC02 are as follows:

Kit batch code	Negative Control (Cat 4)	Negative Control (Cat 20)	Negative Control (Cat 21)
KY25110603	-	-	-
KY25111203	-	-	-
KY25111403	-	-	-

9. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC03 are as follows:

Kit batch code	Negative Control (Cat 14)	Negative Control (Cat 17)	Negative Control (Cat 7)
KY25110707	-	-	-
KY25111213	-	-	-
KY25111413	-	-	-

10. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC04 are as follows:

Kit batch code	Negative Control (Cat14)	Negative Control (Cat17)	Negative Control (Cat7)
KY25111001	-	-	-
KY25111215	-	-	-
KY25111415	-	-	-

11. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC05 are as follows:

Kit batch code	Negative Control (Cat7)	Negative Control (Cat1)	Negative Control (Cat14)
KY25110701	-	-	-
KY25111217	-	-	-
KY25111417	-	-	-

12. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC06 are as follows:

Kit batch code	Negative Control (Cat7)	Negative Control (Cat1)	Negative Control (Cat14)
KY25110703	-	-	-
KY25111301	-	-	-
KY25111501	-	-	-

13. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC07 are as follows:

Kit batch code	Negative Control (Cat5)	Negative Control (Cat9)	Negative Control (Cat13)
KY25110705	-	-	-
KY25111303	-	-	-
KY25111503	-	-	-

14. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-FSC08 are as follows:

Kit batch code	Negative Control (Cat5)	Negative Control (Cat1)	Negative Control (Cat14)
KY25111003	-	-	-
KY25111305	-	-	-
KY25111505	-	-	-

15. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX01 are as follows:

Kit batch code	Positive Control (Cat24)	Positive Control (Cat8)	Positive Control (Cat10)
KY25111005	T5+	T5+	T5+
KY25111307	T5+	T5+	T5+
KY25111507	T5+	T5+	T5+

16. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX02 are as follows:

Kit batch code	Negative Control (Cat5)	Negative Control (Cat9)	Negative Control (Cat13)
----------------	-------------------------	-------------------------	--------------------------

KY25111007	-	-	-
KY25111309	-	-	-
KY25111509	-	-	-

17. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX03 are as follows:

Kit batch code	Negative Control (Cat5)	Negative Control (Cat9)	Negative Control (Cat13)
KY25111101	-	-	-
KY25111311	-	-	-
KY25111511	-	-	-

18. The precision (inter-batch variation) results of the Feline allergen-specific antibody rapid test kit G-CAT-MIX04 are as follows:

Kit batch code	Negative Control (Cat4)	Negative Control (Cat20)	Negative Control (Cat21)
KY25111103	-	-	-
KY25111313	-	-	-
KY25111513	-	-	-

19. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC01 are as follows:

Kit batch code	Positive Control (Dog60)	Positive Control (Dog 113)	Positive Control (Dog 147)
KY25111106	T5+	T2+ T3+ T5+	T3+++ T5+
KY25111316	T5+	T2+ T3+ T5+	T3+++ T5+
KY25111516	T5+	T2+ T3+ T5+	T3+++ T5+

Note: G-DOG-ISC01 T2 *Candida albicans*; T3 *Penicillium chrysogenum*; T5 *D. farinae*

The notation "+++" indicates "strongly positive for the allergens in this panel", the same applies hereafter.

20. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC02 are as follows:

Kit batch code	Positive Control (Dog60)		
KY25110602	T5+	T5+	T5+
KY25111202	T5+	T5+	T5+
KY25111402	T5+	T5+	T5+

21. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC03 are as follows:

Kit batch code	Positive Control (Dog60)		
KY25110606	T5+	T5+	T5+
KY25111206	T5+	T5+	T5+
KY25111406	T5+	T5+	T5+

22. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC04 are as follows:

Kit batch code	Negative Control (Dog158)		
KY25110608	-	-	-
KY25111208	-	-	-
KY25111408	-	-	-

23. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC05 are as follows:

Kit batch code	Negative Control (Dog2)		
KY25110502	-	-	-
KY25111210	-	-	-
KY25111410	-	-	-

24. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-ISC06 are as follows:

Kit batch code	Positive Control (Dog2)		
----------------	-------------------------	--	--

KY25110610	T3+	T3+	T3+
KY25111212	T3+	T3+	T3+
KY25111412	T3+	T3+	T3+

25. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC01 are as follows:

Kit batch code	Negative Control (Dog102)		
KY25111108	-	-	-
KY25111318	-	-	-
KY25111518	-	-	-

26. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC02 are as follows:

Kit batch code	Negative Control (Dog155)		
KY25110604	-	-	-
KY25111204	-	-	-
KY25111404	-	-	-

27. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC03 are as follows:

Kit batch code	Negative Control (Dog155)		
KY25110606	-	-	-
KY25111206	-	-	-
KY25111406	-	-	-

28. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC04 are as follows:

Kit batch code	Negative Control (Dog2)		
KY25111002	-	-	-

KY25111216	-	-	-
KY25111416	-	-	-

29. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC05 are as follows:

Kit batch code	Negative Control (Dog 5)		
KY25110702	-	-	-
KY25111218	-	-	-
KY25111418	-	-	-

30. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC06 are as follows:

Kit batch code	Positive Control (Dog 102)		
KY25110704	T3+	T3+	T3+
KY25111302	T3+	T3+	T3+
KY25111502	T3+	T3+	T3+

31. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC07 are as follows:

Kit batch code	Negative Control (Dog158)		
KY25110706	-	-	-
KY25111304	-	-	-
KY25111504	-	-	-

32. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-FSC08 are as follows:

Kit batch code	Positive Control (Dog158)		
KY25111004	T5+	T5+	T5+
KY25111306	T5+	T5+	T5+

KY25111506	T5+	T5+	T5+
------------	-----	-----	-----

33. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX01 are as follows:

Kit batch code	Positive Control (Dog60)		
KY25111006	T1+T3+T5+	T1+T3+T5+	T1+T3+T5+
KY25111308	T1+T3+T5+	T1+T3+T5+	T1+T3+T5+
KY25111508	T1+T3+T5+	T1+T3+T5+	T1+T3+T5+

34. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX02 are as follows:

Kit batch code	Positive Control (Dog 102)		
KY25111008	T5+	T5+	T5+
KY25111310	T5+	T5+	T5+
KY25111510	T5+	T5+	T5+

35. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX03 are as follows:

Kit batch code	Positive Control (Dog2)		
KY25111102	T5+	T5+	T5+
KY25111312	T5+	T5+	T5+
KY25111512	T5+	T5+	T5+

36. The precision (inter-batch variation) results of the Canine allergen-specific antibody rapid test kit G-DOG-MIX04 are as follows:

Kit batch code	Negative Control (Dog158)		
KY25111104	-	-	-
KY25111314	-	-	-
KY25111514	-	-	-

Based on the above experimental results, the Feline Allergen-Specific Antibody Rapid Test Kit demonstrates consistent detection results across different batches in accordance with the quality control material specifications, and uniform color development is observed among repeated tests of the same quality control material. Similarly, the Canine Allergen-Specific Antibody Rapid Test Kit shows consistent detection results across different batches in accordance with the quality control material specifications, and uniform color development is observed among repeated tests of the same quality control material. Therefore, there is no inter-batch variation, and the kits meet the required specifications.

2.5 Shelf-Life (Accelerated Stability)

2.5.1 Experimental objective

This experiment was conducted to evaluate the shelf life of the company's Canine/Feline allergen-specific antibody rapid test kit.

2.5.2 Experimental materials:

Product name	Catalogue No.	Kit batch codes
Canine allergen-specific antibody rapid test kit	G-DOG-ISC01	KY25111106 , KY25111316 , KY25111516
	G-DOG-ISC02	KY25110602, KY25111202, KY25111402
	G-DOG-ISC03	KY25110606, KY25111206, KY25111406
	G-DOG-ISC04	KY25110608, KY25111208, KY25111408
	G-DOG-ISC05	KY25110502, KY25111210, KY25111410
	G-DOG-ISC06	KY25110610, KY25111212, KY25111412
	G-DOG-FSC01	KY25111108 , KY25111318 , KY25111518
	G-DOG-FSC02	KY25110604, KY25111204, KY25111404
	G-DOG-FSC03	KY25110606, KY25111206, KY25111406
	G-DOG-FSC04	KY25111002, KY25111216, KY25111416
	G-DOG-FSC05	KY25110702, KY25111218, KY25111418
	G-DOG-FSC06	KY25110704, KY25111302, KY25111502
	G-DOG-FSC07	KY25110706, KY25111304, KY25111504
	G-DOG-FSC08	KY25111004, KY25111306, KY25111506
	G-DOG-MIX01	KY25111006, KY25111308, KY25111508
	G-DOG-MIX02	KY25111008, KY25111310, KY25111510
	G-DOG-MIX03	KY25111102, KY25111312, KY25111512
	G-DOG-MIX04	KY25111104, KY25111314, KY25111514

Product name	Catalogue No.	Kit batch codes
Feline allergen-specific antibody rapid test kit	G-CAT-ISC01	KY25111105 , KY25111315 , KY25111515
	G-CAT-ISC02	KY25110601, KY25111201, KY25111401
	G-CAT-ISC03	KY25110605, KY25111205, KY25111405
	G-CAT-ISC04	KY25110607, KY25111207, KY25111407
	G-CAT-ISC05	KY25110501, KY25111209, KY25111409
	G-CAT-ISC06	KY25110609, KY25111211, KY25111411
	G-CAT-FSC01	KY25111107 , KY25111317 , KY25111517
	G-CAT-FSC02	KY25110603, KY25111203, KY25111403
	G-CAT-FSC03	KY25110707, KY25111213, KY25111413
	G-CAT-FSC04	KY25111001, KY25111215, KY25111415
	G-CAT-FSC05	KY25110701, KY25111217, KY25111417
	G-CAT-FSC06	KY25110703, KY25111301, KY25111501
	G-CAT-FSC07	KY25110705, KY25111303, KY25111503
	G-CAT-FSC08	KY25111003, KY25111305, KY25111505
	G-CAT-MIX01	KY25111005, KY25111307, KY25111507
	G-CAT-MIX02	KY25111007, KY25111309, KY25111509

	G-CAT-MIX03	KY25111101, KY25111311, KY25111511
	G-CAT-MIX04	KY25111103, KY25111313, KY25111513

1. Quality control materials

Control ID	Quality control attributes	Properties of quality control materials
Cat 1	Negative Control	<i>Alternaria alternata</i> (-), <i>Malassezia</i> (-), <i>Tyrophagus putrescentiae</i> (-), <i>Blomia tropicalis</i> (-), <i>D. farinae</i> (+)
Cat 4	Negative Control	Flea (-), Mosquito (-), Chicken feather(-), Duck feather(-), Goose feather(-), White ash (-), Mulberry (-), Cottonwood (-), Willow (-), Sycamore (-), Cocklebur (-), Nettle (-), Japanese hop (-), Common ragweed (-), Mugwort (-), Soybean (-), Chicken (-), Duck (-), Sardine (-), Salmon (-), Lamb (-), Pork (-), Turkey meat (-), Shrimp (-), Scallop (-), Squid (-), Wheat (-), Rice (-), Corn (-), Sesame (-), Peanut (-), Cashew (-), Almond (-), Pistachio (-), Walnut (-), Pumpkin (-), Potato (-), Carrot (-), Cauliflower (-), Tomato (-), Apple (-), Pear (-), Pineapple (-), Watermelon (-), Kiwi (-), <i>Blomia tropicalis</i> (-), <i>Malassezia</i> (-), <i>Alternaria alternata</i> (-)
Cat 8	Positive Control	Tree grass pollen (-), Mold (-), Flea (-), Animal dander (-), <i>D. farinae</i> (+)
Cat 9	Positive Control	<i>Alternaria alternata</i> (-), <i>Malassezia</i> (-), <i>Tyrophagus putrescentiae</i> (-), <i>Blomia tropicalis</i> (-), <i>D. farinae</i> (+)
Cat 10	Positive Control	Mugwort (-), <i>Candida albicans</i> (-), <i>Penicillium chrysogenum</i> (-), Cockroach (-), <i>D. farinae</i> (++)
Cat 24	Positive Control	Soybean (-), Sardine (-), Chicken (-), Flea (-), <i>D. farinae</i> (+)

Control ID	Quality control attributes	Properties of quality control materials
Dog2	Positive Control	White ash (-), Mulberry (-), Cottonwood (-), Willow (-), Sycamore (-), Cocklebur (-), Nettle (-), Japanese hop (+), Common ragweed (+), Mugwort (+), Shrimp (-), Scallop (-), Squid (-), Sardine (-), Salmon (-), Turkey meat (-), <i>Malassezia</i> (-)
Dog 5	Positive Control	Shrimp (-), Scallop (+), Squid (-), Sardine (-), Salmon (-)
Dog14	Positive Control	Tree grass pollen (-), Mold (-), Flea (-), Animal dander (-), <i>D. farinae</i> (+), <i>Alternaria alternata</i> (-), <i>Malassezia</i> (-), <i>Tyrophagus putrescentiae</i> (-), <i>Blomia tropicalis</i> (+)
Dog 27	Negative Control	Tree grass pollen (-), Mold (-), Flea (-), Animal dander (-), <i>D. farinae</i> (-), <i>Alternaria alternata</i> (-),

		Malassezia (-), Tyrophagus putrescentiae (-), Blomia tropicalis (-)
Dog 99	Positive Control	Tree grass pollen (-), Mold (-), Flea (-), Animal dander (-), D. farinae (+), Alternaria alternata (-), Malassezia (-), Tyrophagus putrescentiae (-), Blomia tropicalis (-)
Dog 102	Positive Control	Wheat (-), Salmon (-), Duck (-), Mosquito (-), Blomia tropicalis (+)
Dog 155	Negative Control	Flea (-), Mosquito (-), Chicken feather(-), Duck feather(-), Goose feather(-), Squid (-), Lamb (-), Alternaria alternata (-), Mugwort (-), Sycamore (-)
Dog 158	Positive Control	Apple (-), Pear (-), Pineapple (-), Watermelon (-), Kiwi (+)
Dog 162	Negative Control	Lamb (-), Pork (-), Turkey meat (-), Duck (-), Chicken (-), Pumpkin (-), Potato (-), Carrot (-), Cauliflower (-), Tomato (-)
Dog166	Positive Control	Soybean (-), Chicken (-), Duck (-), Sardine (-), Salmon (-), Wheat (-), Rice (-), Corn (-), Soybean (-), Sesame (-), Peanut (-), Cashew (-), Almond (+), Pistachio (-), Walnut (-), Milk (+)

2.5.3 Experimental Method (n=12):

The experiment was conducted according to the instructions. Three batches of the kit were tested. The kits were stored at 37 ° C for 7 days, 14 days, and 20 days, respectively, and tested on the same day. Three replicate tests were performed for each temperature condition. A total of three batches were tested, yielding nine sets of data, to evaluate the shelf life of the kit.

The specific experimental data results are as follows:

1. Feline allergen-specific antibody rapid test kit G-CAT-ISC01; Stability results as follows:

Kit batch code	Positive Control (Cat 8)			Positive Control (Cat 10)			Positive Control (Cat 24)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
KY25111105	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
KY25111315	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+
KY25111515	T5+	T5+	T5+	T5++	T5++	T5++	T5+	T5+	T5+

2. Feline allergen-specific antibody rapid test kit G-CAT-ISC02; Stability results as follows:

Kit batch code	Positive Control (Cat 8)		
	7 days	14 days	20 days
KY25110601	T5+	T5+	T5+
KY25111201	T5+	T5+	T5+
KY25111401	T5+	T5+	T5+

3. Feline allergen-specific antibody rapid test kit G-CAT-ISC03; Stability results as follows:

Kit batch code	Positive Control (Cat 1)			Positive Control (Cat 9)		
	7 days	14 days	20 days	7 days	14 days	20 days
KY25110605	T5+	T5+	T5+	T5+	T5+	T5+
KY25111205	T5+	T5+	T5+	T5+	T5+	T5+
KY25111405	T5+	T5+	T5+	T5+	T5+	T5+

4. Feline allergen-specific antibody rapid test kit G-CAT-ISC04; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110607	-	-	-

KY25111207	-	-	-
KY25111407	-	-	-

5. Feline allergen-specific antibody rapid test kit G-CAT-ISC05; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110501	-	-	-
KY25111209	-	-	-
KY25111409	-	-	-

6. Feline allergen-specific antibody rapid test kit G-CAT-ISC06; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110609	-	-	-
KY25111211	-	-	-
KY25111411	-	-	-

7. Feline allergen-specific antibody rapid test kit G-CAT-FSC01; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25111107	-	-	-
KY25111317	-	-	-
KY25111517	-	-	-

8. Feline allergen-specific antibody rapid test kit G-CAT-FSC02; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110603	-	-	-
KY25111203	-	-	-
KY25111403	-	-	-

9. Feline allergen-specific antibody rapid test kit G-CAT-FSC03; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110707	-	-	-
KY25111213	-	-	-
KY25111413	-	-	-

10. Feline allergen-specific antibody rapid test kit G-CAT-FSC04; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25111001	-	-	-
KY25111215	-	-	-
KY25111415	-	-	-

11. Feline allergen-specific antibody rapid test kit G-CAT-FSC05; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110701	-	-	-
KY25111217	-	-	-

KY25111417	-	-	-
------------	---	---	---

12. Feline allergen-specific antibody rapid test kit G-CAT-FSC06; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110703	-	-	-
KY25111301	-	-	-
KY25111501	-	-	-

13. Feline allergen-specific antibody rapid test kit G-CAT-FSC07; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25110705	-	-	-
KY25111303	-	-	-
KY25111503	-	-	-

14. Feline allergen-specific antibody rapid test kit G-CAT-FSC08; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25111003	-	-	-
KY25111305	-	-	-
KY25111505	-	-	-

15. Feline allergen-specific antibody rapid test kit G-CAT-MIX01; Stability results as follows:

Kit batch code	Positive Control (Cat 24)
----------------	---------------------------

	7 days	14 days	20 days
KY25111005	T5+	T5+	T5+
KY25111307	T5+	T5+	T5+
KY25111507	T5+	T5+	T5+

16. Feline allergen-specific antibody rapid test kit G-CAT-MIX02; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25111007	-	-	-
KY25111309	-	-	-
KY25111509	-	-	-

17. Feline allergen-specific antibody rapid test kit G-CAT-MIX03; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25111101	-	-	-
KY25111311	-	-	-
KY25111511	-	-	-

18. Feline allergen-specific antibody rapid test kit G-CAT-MIX04; Stability results as follows:

Kit batch code	Negative Control (Cat 4)		
	7 days	14 days	20 days
KY25111103	-	-	-
KY25111313	-	-	-
KY25111513	-	-	-

19. Canine allergen-specific antibody rapid test kit G-DOG-ISC01; Stability results as follows:

Kit batch code	Positive Control (Dog 14)			Negative Control (Dog 27)			Positive Control (Dog 99)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
KY25111106	T2+T3+T5++	T2+T3+T5++	T2+T3+T5++	T3+T5+	T3+T5+	T3+T5+	T5++	T5++	T5++
KY25111316	T2+T3+T5++	T2+T3+T5++	T2+T3+T5++	T3+T5+	T3+T5+	T3+T5+	T5++	T5++	T5++
KY25111516	T2+T3+T5++	T2+T3+T5++	T2+T3+T5++	T3+T5+	T3+T5+	T3+T5+	T5++	T5++	T5++

20. Canine allergen-specific antibody rapid test kit G-DOG-ISC02; Stability results as follows:

Kit batch code	Positive Control (Dog 14)			Negative Control (Dog 27)			Positive Control (Dog 99)		
	7 days	14 days	20 days	7 days	14 days	20 days	7 days	14 days	20 days
KY25110602	T5+	T5+	T5+	-	-	-	T5+	T5+	T5+
KY25111202	T5+	T5+	T5+	-	-	-	T5+	T5+	T5+
KY25111402	T5+	T5+	T5+	-	-	-	T5+	T5+	T5+

21. Canine allergen-specific antibody rapid test kit G-DOG-ISC03; Stability results as follows:

Kit batch code	Positive Control (Dog 14)			Negative Control (Dog 27)			Positive Control (Dog 99)		
	7 days	7 days	14 days	20 days	7 days	14 days	20 days	14 days	20 days
KY25110606	T4+T5+	T4+T5+	T4+T5+	-	-	-	T5+	T5+	T5+
KY25111206	T4+T5+	T4+T5+	T4+T5+	-	-	-	T5+	T5+	T5+
KY25111406	T4+T5+	T4+T5+	T4+T5+	-	-	-	T5+	T5+	T5+

22. Canine allergen-specific antibody rapid test kit G-DOG-ISC04; Stability results as follows:

Kit batch code	Negative Control (Dog 155)		
	7 days	14 days	20 days

KY25110608	-	-	-
KY25111208	-	-	-
KY25111408	-	-	-

23. Canine allergen-specific antibody rapid test kit G-DOG-ISC05; Stability results as follows:

Kit batch code	Negative Control (Dog 16)		
	7 days	14 days	20 days
KY25110502	-	-	-
KY25111210	-	-	-
KY25111410	-	-	-

24. Canine allergen-specific antibody rapid test kit G-DOG-ISC06; Stability results as follows:

Kit batch code	Positive Control (Dog 2)		
	7 days	14 days	20 days
KY25110610	T3+T4+T5+	T3+T4+T5+	T3+T4+T5+
KY25111212	T3+T4+T5+	T3+T4+T5+	T3+T4+T5+
KY25111412	T3+T4+T5+	T3+T4+T5+	T3+T4+T5+

25. Canine allergen-specific antibody rapid test kit G-DOG-FSC01; Stability results as follows:

Kit batch code	Positive Control (Dog 166)		
	KY25111108		
	7 days	14 days	20 days

	T1 Soybean	T2 Egg white	T3 Beef	T4 Shrimp	T5 Milk	T1 Soybean	T2 Egg white	T3 Beef	T4 Shrimp	T5 Milk	T1 Soybean	T2 Egg white	T3 Beef	T4 Shrimp	T5 Milk
KY25111108	-	-	-	-	+	-	-	-	-	+	-	-	-	-	+
KY25111318	-	-	-	-	+	-	-	-	-	+	-	-	-	-	+
KY25111518	-	-	-	-	+	-	-	-	-	+	-	-	-	-	+

26. Canine allergen-specific antibody rapid test kit G-DOG-FSC02; Stability results as follows:

Kit batch code	Positive Control (Dog 166)		
	7 days	14 days	20 days
KY25110604	-	-	-
KY25111204	-	-	-
KY25111404	-	-	-

23. Canine allergen-specific antibody rapid test kit G-DOG-FSC03; Stability results as follows:

Kit batch code	Negative Control (Dog 162)		
	7 days	14 days	20 days
KY25110606	-	-	-
KY25111206	-	-	-
KY25111406	-	-	-

24. Canine allergen-specific antibody rapid test kit G-DOG-FSC04; Stability results as follows:

Kit batch code	Positive Control (Dog 5)		
	7 days	14 days	20 days

KY25111002	T2+	T2+	T2+
KY25111216	T2+	T2+	T2+
KY25111416	T2+	T2+	T2+

25. Canine allergen-specific antibody rapid test kit G-DOG-FSC05; Stability results as follows:

Kit batch code	Positive Control (Dog 166)		
	7 days	14 days	20 days
KY25110702	T4+	T4+	T4+
KY25111218	T4+	T4+	T4+
KY25111418	T4+	T4+	T4+

26. Canine allergen-specific antibody rapid test kit G-DOG-FSC06; Stability results as follows:

Kit batch code	Positive Control (Dog 166)		
	7 days	14 days	20 days
KY25110704	T3+	T3+	T3+
KY25111302	T3+	T3+	T3+
KY25111502	T3+	T3+	T3+

27. Canine allergen-specific antibody rapid test kit G-DOG-FSC07; Stability results as follows:

Kit batch code	Negative Control (Dog 162)		
	7 days	14 days	20 days
KY25110706	-	-	-
KY25111304	-	-	-
KY25111504	-	-	-

28. Canine allergen-specific antibody rapid test kit G-DOG-FSC08; Stability results as follows:

Kit batch code	Positive Control (Dog 158)		
	7 days	14 days	20 days
KY25111004	T5+	T5+	T5+
KY25111306	T5+	T5+	T5+
KY25111506	T5+	T5+	T5+

29. Canine allergen-specific antibody rapid test kit G-DOG-MIX01; Stability results as follows:

Kit batch code	Positive Control (Dog 14)		
	7 days	14 days	20 days
KY25111006	T5+	T5+	T5+
KY25111308	T5+	T5+	T5+
KY25111508	T5+	T5+	T5+

30. Canine allergen-specific antibody rapid test kit G-DOG-MIX02; Stability results as follows:

Kit batch code	Positive Control (Dog 102)		
	7 days	14 days	20 days
KY25111008	T5+	T5+	T5+
KY25111310	T5+	T5+	T5+
KY25111510	T5+	T5+	T5+

31. Canine allergen-specific antibody rapid test kit G-DOG-MIX03; Stability results as follows:

Kit batch code	Positive Control (Dog 2)		
	7 days	14 days	20 days
KY25111102	T4+	T4+	T4+

KY25111312	T4+	T4+	T4+
KY25111512	T4+	T4+	T4+

32. Canine allergen-specific antibody rapid test kit G-DOG-MIX04; Stability results as follows:

Kit batch code	Negative Control (Dog 155)		
	7 days	14 days	20 days
KY25111104	-	-	-
KY25111314	-	-	-
KY25111514	-	-	-

The experimental results demonstrate that the test outcomes meet the requirements for both positive and negative controls, with consistent color development across repeated tests of the same control sample. In accordance with the product technical specifications, the kit must be stored at 37°C for 20 days from the date of production (with a shelf life of 24 months). Consequently, The following kits have a shelf life of 24 months: Feline allergen-specific antibody rapid test kit , Canine allergen-specific antibody rapid test kit.